# IDROSSILAPATITE AT<sup>®</sup> CE

0373

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# IDROSSILAPATITE AT<sup>®</sup> CE

0373

Granular osteoconductive Hydroxylapatite for guided tissue regeneration in dental surgery

"Guided regeneration materials system  
Collagene AT<sup>®</sup> Idrossilapatite AT<sup>®</sup>

## Properties and Technical Features

The **Idrossilapatite AT<sup>®</sup>** application prevents epithelial and connective tissue cells from interfering with bone growth and periodontal tissue growth.

The maintenance of a space beneath the membrane by **Idrossilapatite AT<sup>®</sup>** allows production of new bone tissue and new periodontal tissue and for the prevention of post-surgical reabsorption of the alveolar ridges.

Surface nanometric micro-irregularity of **Idrossilapatite AT<sup>®</sup>** allows the absorption of plasmatic proteins with function of growth factors and with osteoinductive action.

**Idrossilapatite AT<sup>®</sup>**, calcium hydroxylapatite with molecular weight 502,32 with nanometric microgranulometry at surface, has been designed and created because its particular granulometry allows the direct manipulation and application of the material into the operating field, thanks to its reciprocal adhesion of the particles.

Sphere shaped granulometry 200-500 micron adhere each other for capillarity and leave continuous intermediate spaces, which have an osteoconductive action and are progressively colonized by bony trabeculations.

Reabsorbable progressively in centripetal direction, towards the inside of the granule in a time span of 6-60 months, during the formation of new bone.

The absorption time corresponds to the physiological time for the trajectorial remodeling of trabecula bone.

To get a good result in Guided Tissue Regeneration, it must be placed **Idrossilapatite AT<sup>®</sup>** under the Collagene AT<sup>®</sup> membrane.

## Packaging and Conservation

**Idrossilapatite AT<sup>®</sup>** are sold in three different packagings:

- 1 glass bottle of 17 g with bake-lite cap. **(code 2600).**
- 10 glass bottle single use of 1 g with cap **(code 2601).**
- 2 glass bottle single use of 1 g with screw cap in welded blister **(code 2609).**

**Idrossilapatite AT<sup>®</sup>** is sterilized by gamma radiation.

## Main Utilizations

### In Implant Surgery

allows the application of regeneration materials to prevent post-operative ridge atrophy and atrophic ridges augmentation before implant surgery.

### In Periodontal Surgery

allows guided regeneration in intrabony pockets.

### In Endodontic Surgery

allows new bone growth in the periapical zone without residual cavity.

### In Oral Surgery

allows the filling of post-extractive sockets and prevents bone atrophy.

## Precautions

**Idrossilapatite AT<sup>®</sup>** is contraindicated in patients with acute infections or contaminated wounds.

Medical device **Idrossilapatite AT<sup>®</sup>** can be used only by Dental Surgeons, Physicians, Dentists qualified in surgical dental procedures.

Before using **Idrossilapatite AT<sup>®</sup>** make sure the packaging is intact.

The including adhesive tracking labels shall appear in the patient data.

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