



CENTRO DI ODONTOIATRIA OPERATIVA s.r.l.
Divisione Materiali: Via Guizza, 309 - 35125 Padova
Tel. 049.681123 - Fax 049.8806560 - P.IVA 01385450281
<http://www.odopguizza.it> - e-mail: info@odopguizza.it

TECHNICAL SHEET

Product name: **COLLAGENE AT®**

MEDICAL DEVICE COMPONENT:
"SYSTEM OF MATERIALS FOR GUIDED REGENERATION
COLLAGENE AT® IDROSSILAPATITE AT®"



- Material** > Lyophilized collagen of equine extraction.
- Structure** > Membrane with nanometric micro-irregularities surface and about 20 micron wide porosity, allowing spontaneous adhesion on tissue without sutures and barrier effect against cellular migration.
- Reabsorption** > Within 180 days.
- Sterilization** > Gamma rays.
- Package contents** > 6 quadrangular membranes of 22 x 22 mm (code 4101) in individually sterilized envelopes.
> 6 circular membranes of 17 mm (code 4102) with individually sterilized envelopes.
This form naturally adapts to the convexity of teeth in intraosseous periodontal pockets.
> 3 quadrangular membranes of 22 x 22 mm and 3 circular membranes of 17 mm (code 4103) in individually sterilized envelopes.
- Indications** > Manipulate the membrane with perfectly dry instruments and press with dry gauzes.
The membrane should be applied on Hydroxylapatite AT® used as an osteo-conductor support material or on bone graft to create a perfect tent barrier effect.
- Warnings** > Collagene AT® is single use only.
Supplied in heat-sealed envelopes.
Open only before use.
Store at room temperature.
Avoid exposure to temperatures over 40°.
Use only intact envelopes.
- Packaging** > Cardboard boxes.
- Main Uses** > Guided tissue regeneration in periodontology.
> Guided bone tissue regeneration in implantology.
> Prevents post-extractive atrophy in exodontial.

CENTRO DI ODONTOIATRIA OPERATIVA